



INNOVATIVE MULTI-LUMEN TUBING FOR MINIMALLY INVASIVE AND CARDIAC APPLICATIONS

HelixTwist™ is a proprietary multi-lumen silicone tubing technology designed to prevent kinks or blocked lumens by exerting an equal balance of stress across the inner and outer lumens as a tube bends. HelixTwist™ tubing is optimal for pacemakers, breathing tubes, neurostimulation and other applications that require navigation through tortuous pathways within the human body.

ABOUT FREUDENBERG MEDICAL

Freudenberg Medical is an experienced global partner for the design, development, and manufacture of innovative medical device technologies. Our deep experience and capabilities range from minimally invasive, catheter and handheld technologies to the development and production of challenging silicone, thermoplastic and metal components. We process advanced materials with cutting-edge processes such as precision molding, complex extrusion, laser machining and custom coatings.

Freudenberg Medical Europe GmbH
Liebigstr. 2 – 8
67661 Kaiserslautern, Germany

Published by:
Freudenberg Medical
Marketing
Phone: +49 631 5341 7520
info@freudenbergmedical.de
www.freudenbergmedical.de

August 2021 | Item No. FMED026EN0421

- Maximum of 6 inner lumens
- OD up to 6,00 mm +/- 0,10 mm (~ 0.24")
- Available with implant-grade silicone material and /or radiopaque properties
- Durometer: 60 Shore A – 80 Shore A
- Color: Translucent or colored according to customer specification
- Class 8 cleanroom production, ISO 13485 certified
- Packaging: Double PE bag, sealed, label on inner bag. Tubing coils, on spools or cut-to-length according to customer specification
- In-line measurement control of inner and outer diameter with state-of-the-art Helix iMC™ technology
- Project management and feasibility analysis with help of our in-house Application Engineering team
- Validation: Standard validation protocol readily available or according to customized validation requirements
- In-house test labs to support testing of physical, chemical, rheological properties, or with failure analysis